RECOVERY & MOVEMENT ENHANCEMENT TECHNOLOGY™

As far back as I can remember, I have always had a passion for sports and science. From my own experience as an athlete, I studied the science of the body and learned that maintaining a healthy body is essential to an athlete's performance. Several years ago, I set out to create a brand whose products would improve athletic performance by enhancing athletic recovery and preventing injury. In order to develop products of this caliber, I collaborated with some of the world's best athletes and trainers so that the products would meet their standards. Together our goal was to create a series of products that would be a set of tools athletes could use to recover faster from training and competition, and prevent injury by keeping the body in proper balance. I would go on to name the brand after our first product- the Hyperice: the most advanced portable ice compression device on the market. Today our products are used by over one thousand of the top athletes in the world. We strive to perfect all of our products and look forward to collaborating with more professional athletes to develop the world's best prevention and recovery devices on the planet.

Anthony Katz | Founder  
Hyperice, Inc

## CRYOTHERAPY FOR INJURY REPAIR & PERFORMANCE ENHANCEMENT

The use of ice, or cryotherapy, for musculoskeletal injuries has been a primary treatment approach used by healthcare providers for many years. Today, cryotherapy use continues to be a proven and beneficial adjunct to managing soft tissue injury. It is proposed that applying cooling to a musculoskeletal injury decreases the tissue temperature, resulting in diminished pain, cellular metabolism, and muscle spasm, thus minimizing the inflammatory response and improving recovery after soft tissue trauma.

### THE INFLAMMATORY RESPONSE

An inflammatory response occurs at the site of the injury and is the local, tissue-level response to the injury and the mechanism by which a tissue heals. Signs of tissue inflammation may include pain, redness, warmth, swelling, and decreased function. Inflammation can develop from both acute injuries, such as ligament sprains or tears, and chronic injuries, such as arthritis or tendinitis. Swelling of the injured region is not uncommon and contributes to loss of movement and ultimately a decrease in function.

### THERAPEUTIC EFFECTS OF ICE

Cryotherapy, or cold therapy, is a common modality used in the treatment of soft tissue injuries. Cold therapy has an impact on a tissue's metabolism, pain response, and hemodynamics, thus aiding in the management of the inflammatory response. It has also been demonstrated that the application of cryotherapy decreases blood flow to the injured region, thus reducing edema formation. The combined effects of decreased energy requirements, reduced blood flow, analgesic response, and swelling reduction all diminish the tissue's inflammatory response, ultimately resulting in a shortened healing process and hastened return to activity.

### ROLE OF COMPRESSION

Compression is often used in conjunction with cryotherapy. Benefits of compression include improved contact between the skin and the cold source, greater reduction of blood flow to the region, and an increased insulation effect, which may further reduce tissue temperatures. Compression also assists with control of edema formation that may arise after injury or secondary to microtrauma sustained during a hard workout.

EACH HYPERICE DEVICE IS COMPOSED OF TWO PARTS THAT WORK TOGETHER AS A SYSTEM TO MAXIMIZE COLD + COMPRESSION AND DELIVER OPTIMAL COLD THERAPY.

## 1) ICE CELL W/ AIR RELEASE TECHNOLOGY

The Ice Cell provides the “cold” component by harnessing the power of real ice, or HYPERICE Fuel, which is loaded into the Ice Cell, used in conjunction with the patented Air Release Valve, which removes air from the Ice Cell to enable maximum compression throughout the icing process.

### PRECISION HARDWARE

Wide-mouth opening for easy filling and flange collar to stabilize Ice Cell during treatments.

### AIR RELEASE VALVE


Patented Air Release Valve removes air from the Ice Cell for optimal compression.

### ULTRA-THIN SKIN

Ultra-thin antimicrobial skin harnesses the power of real ice to quickly transfer cold to the body.

## 2) COMPRESSION WRAP

Each Compression Wrap is designed to apply firm pressure to the Ice Cell, pushing all excess air trapped in the cell to the surface. By pressing the Air Release Valve, the excess air is removed from the Ice Cell, allowing the Compression Wrap to make direct contact with the ice, forming an ice “cast” around the specific body part.

\*SHOULDER-RIGHT DEVICE SHOWN

## SHOULDER DEVICE

Right Shoulder / Left Shoulder (Sold Separate)

Sleek, compact, and comfortable. For the athlete that throws, swings or strikes. Delivers optimal cold to the shoulder to best experience the benefits provided by cold therapy.

## UTILITY DEVICE

Ankle / Elbow / Shin / Wrist / Youth Knee

The fixer. Designed to go almost anywhere. Compact design fitted with smaller ice cell, for targeted cold and compression. Delivers optimal cold to the ankle, elbow, shin or wrist. Also fits excellent for the knee of youth athletes.

## KNEE DEVICE

Right Knee / Left Knee

Recharge + Repair + Recover. Essential gear for any athlete who moves fast and needs to recover faster. Delivers optimal cold to the knee to best experience the benefits provided by cold therapy.

## BACK DEVICE

Lower Back / Middle Back

Revitalize your core. Restabilize your body. Multi-layered compression for the back. Rebuild from the center, return to balance. Delivers optimal cold to the back to best experience the benefits provided by cold therapy.

VYPER 2.0

## V Y P E R 2 . 0

HIGH-INTENSITY VIBRATING FITNESS ROLLER

The Vyper 2.0 is the next generation of the world's most powerful vibrating fitness roller. Now stronger, sleeker, and even more durable.

To purchase or learn more visit:  
[hyperice.com](http://hyperice.com)

HYPERSPHERE

## HYPERSPHERE

HIGH-INTENSITY VIBRATING FITNESS BALL

Powerful yet small, the Hypersphere targets your tightest areas to unlock the body and help you move better.

To purchase or learn more visit:  
[hyperice.com](http://hyperice.com)

VENOM

HYPERICE FUEL

VENOM

---

WEARABLE HEAT + VIBRATION BACK DEVICE

---

The Venom is a cutting edge, digitally connected wearable device that combines heat and vibration to warmup, loosen and relax muscles.

To purchase or learn more visit:  
[hyperice.com](http://hyperice.com)

HYPERVOLT

---

VIBRATION MASSAGE DEVICE

---

The Hypervolt is a cordless state-of-the-art vibration massage device that helps relax sore and stiff muscles to improve mobility.

To purchase or learn more visit:  
[hyperice.com](http://hyperice.com)

# hypersoothe™

ESSENTIAL OIL-INFUSED SOOTHING CREAM

Hypersoothe is a lightweight topical cream infused with essential oils and natural plant extracts to soothe the body's muscles and joints. Formulated to meet the standards of the world's best athletes.

To purchase or learn more visit:  
[hyperice.com](http://hyperice.com)

# HYPERICE FUEL™

REUSABLE SYNTHETIC ICE

A reusable synthetic ice, Hyperice Fuel™ provides optimal freezer temperatures through a full icing cycle. Hyperice Fuel can be used with any of the Hyperice Ice Compression Devices.

To purchase or learn more visit:  
[hyperice.com](http://hyperice.com)

## INSTRUCTIONS

### USAGE INSTRUCTIONS

**1** Fill the Ice Cell 2/3 with ice – do not add water. Crushed ice or small ice nuggets work best with all HYPERICE® products, however ice cubes may be substituted as needed.

**2** Place the Ice Cell on a flat surface and flatten out - removing as much air from the Ice Cell as possible before securing the cap.

**3** Insert the Ice Cell into the Compression Wrap.

**4** Place on the body and apply the compression straps to the desired level of compression.

**5** Lightly press the Air Release Valve to remove air from the Ice Cell. For added compression, re-tighten the compression straps after the air has been released.

**6** Leave the device on for 15 – 20 minutes. Press the Air Release Valve once every 5 minutes to release any air that accumulates during the treatment to sustain optimal compression.

### CARE INSTRUCTIONS

**1** After use, rinse the Ice Cell inside and out with fresh water and air dry.

**2** After the Ice Cell is dry, store in a cool, dry place.

**3** Do not expose the Ice Cell to direct sunlight and/or heat, as the UV rays can damage the medical-grade material.

**4** Do not fill the Ice Cell with anything other than ice or HYPERICE FUEL.

**5** Do not machine wash or iron.

### CAUTION

**1** THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS.

**2** DO NOT LEAVE ANY HYPERICE® DEVICE ON THE BODY FOR MORE THAN 25 MINUTES AT A TIME.

For more detailed directions on how to use your HYPERICE® device, please visit the PRODUCTS section of our website at [www.hyperice.com](http://www.hyperice.com). There you will find an instructional video that clearly illustrates how to use your HYPERICE® device.

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## 1-YEAR WARRANTY

If within one year from the date of purchase, this product fails due to a defect in material or workmanship, Hyperice Inc., and/or the international distributor, dealer or retailer from whom the purchase was made, will repair or replace the product, or necessary components, free of charge.

This warranty excludes:

- (a) damage caused by accident, abuse, mishandling, or transport;
- (b) units subjected to unauthorized repair;
- (c) units not used in accordance with Hyperice care instructions;
- (d) damage exceeding the cost of the product;
- (e) deterioration of the delivered product resulting from abnormal storage and/or safeguarding conditions on the client's premises, and
- (f) failure to provide the dated proof of purchase.

Some states and countries do not allow a limitation of damages, so the foregoing limitation may not apply to you. This warranty guarantees specific legal rights, and but other rights may vary from country to country and from state to state. This warranty service is available through the operating policies and procedures of Hyperice Inc., and all Hyperice international distributors.

If you are a client, please contact the appropriate Hyperice business entity, dealer or retailer, from whom you directly purchased the product for warranty and return authorization procedures. This warranty is initiated and executed by Hyperice Inc.

15440 Laguna Canyon Road, Suite 230  
Irvine, California 92618  
**HYPERICE.COM**